



Fast, Flexible,

Functional & Affordable

Introducing the all-new FlexiStar standard medical cart, offering you out-of-the-box functionality, faster turnaround, smaller minimum quantities and an affordable price point without compromising quality or performance.



### FlexiStar Optional Features

Enjoy features and functionalities often reserved for custom applications without the custom price, large quantities, or lead time.

- Built to ISO 60601 standards
- 3 adjustable shelves
- Spacious bin
- Ergonomic tube handles
- 4 dual 5" locking casters (w/optional cable pushers)
- Cable management
- Keyboard tray
- Cord wraps
- Removable back panel for quick access without tools
- 1-2 drawers
- Single or dual monitor mount (swivel, tilt, pan, and height adjustable)
- Probe holder
- IV pole



10 East Park Ave, Kiel WI 53042

920-894-2233

[www.medicalcarts.org](http://www.medicalcarts.org)

## Customized Options

FlexiStar models feature many standard features to meet most needs, yet can be configured to meet your exacting specifications, including custom:

- Mounting holes for devices
- Variety of articulating arms
- Cylinder tank holder brackets
- Foot switch
- Machined or casted metal features
- Plastic features including:
  - Thermoforming
  - Injection molding
  - Cast urethane
  - 3D printing
- And More...



## FlexiStar Baseline Dimensions

- Height 39"
- Base: 21.75" x 28.5"
- Lower shelf: 13.25" x 20.25"
- Middle shelf: 18.5" x 20.5"
- Top surface: 20.5" x 20.5"
- Monitor size: 19"-24"
- Height adjustable monitor: 6"
- Cart weight: 100lbs
- Casters: 5"

Please Contact  
HUI's Medical Business  
Development Team  
1-800-877-8913  
MBD@huimfg.com

