

VERSAPLUS™

Brought to you by HUI

Point-of-care Powerhouse.

Looking for a workstation that really works? Configure your cart to create the right workspace for the job.

A robust workspace in a practical package.

Our VersaPlus™ standard cart offers plenty of workspace in a versatile cart. Our rear-column design optimizes space for surface monitors or tablets, with basket, shelf, or lockable drawer options below so staff can manage tools and accessories in an organized system.

Cables can be managed inside the column or with installed cable wraps. Choose a large or small top surface options to fit your device monitor or tablet.

VersaPlus™ is designed to maximize modular configurations while not compromising on safety and ergonomics. The four 5" medical-grade casters maximize freedom of movement and stability at the point of care.



Quality Construction

Built to IEC 60601-1 Standards

Manufactured on-site
in our midwest facility

Machined or casted
metal components

Compact Profile

Cart weight under 50 lbs.

24 inch base ideal for
small-scale environments

Rear column cable
management

Diverse Applications

Telemedicine

Diagnostic/Monitoring
Equipment


EEG/Cardiology

ENT Treatment

VERSAPLUS™

Expedite your speed to market.

Contact sales@huimfg today to get started!



VersaPlus Baseline Dimensions:

- Height: 42.7"
- Base: 19.3" x 15.7"
- Cart Weight: 48 lbs.
- Casters: 4 x 5" locking casters

Customizable Features:

- Large or small top surface options for device fit
- Drawer or shelf height and orientation on column
- Cable management in column or additional cord wrap
- Professionally-matched powder-coat color and finish

Accessory Options:

- Add monitor or tablet mount to top work surface
- Small or large basket options
- Option for drawer mounted beneath work surface

Shipping and Assembly:

- Can ship fully-assembled or in flatpacks
- Flatpack option kitted for easy assembly
- Packaging is engineered to minimize damage in transit

